ZENOTEC Color Zr

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English

ZENOTEC Color Zr

ZENOTEC Color Zr is a range of ready-to-use colouring solutions for staining dental frameworks made from zirconium oxide prior to sintering. The shades of ZENOTEC Color Zr are coordinated with ZENOTEC Zr_{Bridge} (white) and ZENOTEC Zr_{Crown+} (white).

ZENOTEC Color Zr_{Fix}

ZENOTEC Color Zr_{Fix} is a range of ready-to-use solutions for fixing the colour of zirconium oxide frameworks shaded with ZENOTEC Color Zr.

Safety first

Avoid direct contact with the skin. In the event of eye contact, rinse open eye well with plenty of water and seek medical advice.

It is recommended to wear gloves, goggles and suitable protective clothing. ZENOTEC Color Zr can stain clothes and sensitive surfaces. Ensure adequate ventilation at the workplace when working with ZENOTEC Color Zr_{Fix}. Always seal containers firmly after use.

How to use

ZENOTEC Color Zr is a range of ready-to-use colouring solutions for staining zirconium dioxide frameworks to match the colours of the standard shade guide. The colouring solutions are used at room temperature. In order to stain the unsintered frameworks, the solutions are poured into the plastic beakers (with basket) supplied. (If any other vessel is used, ensure that this is made of plastic or glass and under no circumstances of metal.) The ZENOTEC Color Zr and ZENOTEC Color Zr_{Fix} liquids must not be mixed together.

Use the baskets supplied and the plastic tweezers (do not use metal tweezers) to carefully immerse the unsintered frameworks in the colouring solution until the workpiece is fully covered with liquid. The basket makes it easy to position the framework in the liquid and to remove it after staining. It is sufficient to immerse the frameworks in the solution very briefly (approx. 10 seconds).

After staining, use the plastic tweezers to remove the frameworks from the basket and use absorbent paper to carefully blot the excess liquid from the frameworks. Then immerse the stained frameworks in the ZENOTEC Color Zr_{Fix} fixing solution (using the plastic tweezers) and leave the work in the solution for two minutes. After removing the frameworks from the fixing solution, dab once again with absorbent paper to carefully blot off the excess liquid.

The stained frameworks **must be completely dry** before they are sintered, otherwise damage to both the frameworks and the sinter furnace may result. As a guideline, the following table gives typical values for drying temperature and drying time:

Drying temperature	80 °C
Drying time, minimum	1 – 3 units: 30 min 4 – 6 units: 2 hours > 6 units: 4 hours

Solid frameworks (e.g. with gingival mask) and orally connected frameworks require a drying time of > 4 hours. The frameworks can be dried in the furnace once it is preheated to drying temperature.

Veneering

Frameworks coloured with ZENOTEC Color Zr are shaded to match the ZIROX veneers. The required tooth shade is attained after veneering. In order to obtain greater overall translucence, no liner should be used.

Storage

- Store ZENOTEC Color Zr and ZENOTEC Color Zr_{Fix} at room temperature in tightly closed containers
- Protect the solutions from frost
- Do not expose to direct sunlight
- Keep solutions as clean as possible
- If the solutions should become cloudy, refrain from further use
- Do not decant or store in metal containers, even for a short time

Waste disposal

The used solutions can be poured down the drain after being diluted with plenty of water. Please observe local waste water disposal regulations.

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GmbH & Co. KG