

Safety Data Sheet

according to 1907/2006/EC, Article 31

Printing date 01.04.2010

Revision: 01.04.2010

1 Identification of the substance/preparation and of the company/undertaking

- **Product details**

- **Trade name: Porta Lot 1090 W**

- **Application of the substance / the preparation**

Precious metall alloy for dental technology

Hard solder

- **Manufacturer/Supplier:**

Wieland Dental + Technik GmbH & Co.KG

Schwenninger Str. 13

75179 Pforzheim

Telefon + 49 (07231)-37050, Telefax + 49 (07231)-357959

- **Further information obtainable from:**

Abteilung Umweltschutz

Marcus.Mueller@wieland-dental.de

- **Information in case of emergency:**

+49 (07231)-37050

(Be kindly advised that this number can only be contacted at usual office hours.)

2 Hazards identification

- **Hazard description:** Not applicable.

- **Information concerning particular hazards for human and environment:**

The product does not have to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

- **Classification system:**

The classification is according to the latest editions of the EU-lists, and extended by company and literature data.

- **GHS label elements**

Safety data sheet available on request.




- **Additional information:** In compact form is this alloy not dangerous for human or the aquatic environment.

3 Composition/information on ingredients

- **Chemical characterization**

- **Description:** Metal in compact form.

- **Dangerous components:**

CAS: 7440-22-4 EINECS: 231-131-3	silver		25-50%
CAS: 7440-05-3 EINECS: 231-115-6	palladium	 O; R 8 Danger:  2.14/2 Warning:  2.11/2	10-25%

- **Additional information:** For the wording of the listed risk phrases refer to section 16.

4 First-aid measures

- **After inhalation:** Supply fresh air; consult doctor in case of complaints.

- **After skin contact:**

After contact with the molten product, cool rapidly with cold water.

If skin irritation continues, consult a doctor.

- **After eye contact:**

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

- **After swallowing:** Seek immediate medical advice.

GB

(Contd. on page 2)

Safety Data Sheet

according to 1907/2006/EC, Article 31

Printing date 01.04.2010

Revision: 01.04.2010

Trade name: Porta Lot 1090 W

(Contd. of page 1)

5 Fire-fighting measures

- **Suitable extinguishing agents:**
CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Protective equipment:** No special measures required.

6 Accidental release measures

- **Person-related safety precautions:** Not required.
- **Measures for environmental protection:** No special measures required.
- **Measures for cleaning/collecting:**
Pick up mechanically.
Dispose contaminated material as waste according to item 13.

7 Handling and storage

- **Handling:**
- **Information for safe handling:**
Extractors are required on all machines used for thermal processing or splinter removal processes.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.

8 Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.

- **Ingredients with limit values that require monitoring at the workplace:**

7440-22-4 silver

WEL | Long-term value: 0.1 mg/m³

- **Additional information:** The lists valid during the making were used as basis.

- **Personal protective equipment:**
- **General protective and hygienic measures:** Wash hands before breaks and at the end of work.
- **Respiratory protection:** Not required.
- **Protection of hands:**



Protective gloves

- **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

- **Penetration time of glove material**

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **Eye protection:** Safety glasses

GB

(Contd. on page 3)

Safety Data Sheet

according to 1907/2006/EC, Article 31

Printing date 01.04.2010

Revision: 01.04.2010

Trade name: Porta Lot 1090 W

(Contd. of page 2)

9 Physical and chemical properties

· General Information

Form:	Solid
Colour:	White
Odour:	Odourless

· Change in condition

Melting point/Melting range: 980 - 1090°C
Boiling point/Boiling range: Undetermined.

· **Flash point:** Not applicable.

· **Self-igniting:** Product is not selfigniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Density:** -

· **Solubility in / Miscibility with water:** Insoluble.

· **Solvent content:**

· **Solids content:** 100.0 %

10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Dangerous reactions** Reacts with strong acids.
- **Dangerous decomposition products:** Toxic metal oxide smoke

11 Toxicological information

- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.

12 Ecological information

- **General notes:** Generally not hazardous for water

13 Disposal considerations

- **Product:**
- **Recommendation** Contact manufacturer for recycling information.

· European waste catalogue

10 00 00	WASTES FROM THERMAL PROCESSES
10 07 00	wastes from silver, gold and platinum thermal metallurgy
10 07 04	other particulates and dust

- **Uncleaned packaging:**
- **Recommendation:** Packaging may be reused or recycled after cleaning.

(Contd. on page 4)

Safety Data Sheet

according to 1907/2006/EC, Article 31

Printing date 01.04.2010

Revision: 01.04.2010

Trade name: Porta Lot 1090 W

(Contd. of page 3)

- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

14 Transport information

- **Maritime transport IMDG:**
- **Marine pollutant:** No

15 Regulatory information

- **Labelling according to EU guidelines:**
Although this product is not subject to identification regulations, we recommend that the safety suggestions are observed.
- **Special labelling of certain preparations:**
Safety data sheet available for professional user on request.
- **National regulations:**
- **Waterhazard class:** Generally not hazardous for water.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant R-phrases**
8 Contact with combustible material may cause fire.
- **Department issuing MSDS:** Abteilung Umweltschutz
- **Contact:** Herr Marcus Müller
- **Abbreviations and acronyms:**
IMDG: International Maritime Code for Dangerous Goods
GHS: Globally Harmonized System of Classification and Labelling of Chemicals

GB