

Material Safety Data Sheet

According to 91/155 EEC

Printing date 13.06.2005

Reviewed on 13.06.2005

*

1 Identification of the substance/preparation and of the company/undertaking

- **Product details**
- **Trade name:** Imagine h. e. Colours-Liquid
- **Application of the substance / the preparation** Auxiliary for dental technology
- **Manufacturer/Supplier:**
Wieland Dental + Technik GmbH & Co.KG
Schwenninger Str. 13
75179 Pforzheim
Telefon + 49 (07231)-37050, Telefax + 49 (07231)-357959
- **Further information obtainable from:** AT (Tel.:+49(0)6007-9176-222

2 Composition/information on ingredients

- **Chemical characterization**
- **Description:** Mixture of substances listed below with nonhazardous additions.
- **Dangerous components:** Void

3 Hazards identification

- **Hazard description:** Not applicable.
- **Information concerning particular hazards for human and environment:** Not applicable.

4 First-aid measures

- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Rinse with warm water.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** Seek medical treatment.

*

5 Fire-fighting measures

- **Suitable extinguishing agents:**
CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Protective equipment:** Wear self-contained respiratory protective device.

6 Accidental release measures

- **Person-related safety precautions:** Not required.
- **Measures for environmental protection:** No special measures required.
- **Measures for cleaning/collecting:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Additional information:** No dangerous substances are released.

*

7 Handling and storage

- **Handling:**
- **Information for safe handling:** Ensure good ventilation/exhaustion at the workplace.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.

(Contd. on page 2)

GB

Material Safety Data Sheet

According to 91/155 EEC

Printing date 13.06.2005

Reviewed on 13.06.2005

Trade name: **Imagine h. e. Colours-Liquid**

(Contd. of page 1)

- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.

*

8 Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.
- **Ingredients with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists valid during the making were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:** Not necessary if room is well-ventilated.
- **Protection of hands:** Not required.
- **Material of gloves:** Void
- **Penetration time of glove material:** Void
- **Eye protection:** Not required.

*

9 Physical and chemical properties

· General Information

Form:	Fluid
Colour:	Colourless
Odour:	Product specific

· Change in condition

Melting point/Melting range:	-50°C
Boiling point/Boiling range:	206°C

· **Flash point:** 121°C

· **Ignition temperature:** 440°C

· **Danger of explosion:** Product does not present an explosion hazard.

· **Density at 20°C:** 1 g/cm³

· Solubility in / Miscibility with water:

Fully miscible.

· **pH-value:** Not applicable

10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:**
No decomposition if used according to specifications.
- **Dangerous reactions:** No dangerous reactions known.
- **Dangerous decomposition products:** No dangerous decomposition products known.

GB

(Contd. on page 3)

Material Safety Data Sheet

According to 91/155 EEC

Printing date 13.06.2005

Reviewed on 13.06.2005

Trade name: **Imagine h. e. Colours-Liquid**

(Contd. of page 2)

11 Toxicological information

- **Acute toxicity:**

- **LD/LC50 values relevant for classification:**

Oral	LD50	18610 mg/kg (Ratte)
------	------	---------------------

- **Primary irritant effect:**

- **on the skin:** No irritant effect.

- **on the eye:** No irritating effect.

- **Sensitization:** No sensitizing effects known.

- **Additional toxicological information:**

When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

The substance is not subject to classification according to the latest version of the EU lists.

12 Ecological information

- **General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water

13 Disposal considerations

- **Product:**

- **Recommendation** Must be specially treated adhering to official regulations.

- **European waste catalogue**

07 00 00	WASTES FROM ORGANIC CHEMICAL PROCESSES
----------	--

07 01 00	wastes from the manufacture, formulation, supply and use (MFSU) of basic organic chemicals
----------	--

07 01 04	other organic solvents, washing liquids and mother liquors
----------	--

- **Uncleaned packaging:**

- **Recommendation:**

Packagings that may not be cleansed are to be disposed of in the same manner as the product.

- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

14 Transport information

- **Land transport ADR/RID (cross-border)**

- **ADR/RID class:** -

- **Maritime transport IMDG:**

- **IMDG Class:** -

- **Marine pollutant:** No

- **Air transport ICAO-TI and IATA-DGR:**

- **ICAO/IATA Class:** -

15 Regulatory information

- **Labelling according to EU guidelines:**

The product is not subject to identification regulations under EU Directives and the Ordinance on Hazardous Materials (German GefStoffV).

Observe the general safety regulations when handling chemicals.

(Contd. on page 4)

Material Safety Data Sheet

According to 91/155 EEC

Printing date 13.06.2005

Reviewed on 13.06.2005

Trade name: Imagine h. e. Colours-Liquid

(Contd. of page 3)

- **National regulations:**
- **Waterhazard class:** Water hazard class 1 (Self-assessment): slightly hazardous for water.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Department issuing MSDS:** Abteilung Umweltschutz
- **Contact:** Herr Selzer
- *** Data compared to the previous version altered.**

GB